

Name: _____ Date: _____

Combine with a Backboard



Combine with a Backboard Practical

DIRECTIONS: Review the Combine with a Backboard video in HILT U. Mark whether the candidate completed the task correctly or incorrectly.

Task	Correct	Incorrect	Notes
<u>Two person lift from ends of backboard</u>			
Laid the backboard on top of a full length HILT.			
Threaded the four low lift base handles into the four backboard corner handholds.			
Packaged the patient according to local directives.			
Used correct lifting technique to minimize injury exposure to self.			
Lifted the patient in unison with team lift.			
Lowered the patient safely.			
<u>Two person lift from sides of backboard</u>			
Folded each end of the HILT up to the second sets of extenders.			
Laid the backboard on top of the shortened HILT. Avoided an imbalance by not combining too close to the backboard head or foot.			
Deployed the second sets of extenders at each corner of the shortened HILT.			
Threaded each extender through the nearest backboard handhold.			
Packaged the patient according to local directives.			
Used correct lifting technique to minimize injury exposure to self.			
Lifted the patient in unison with team lift.			
Lowered the patient safely.			

Name: _____ Date: _____

Combine with a Backboard

Multi-Person Lift with full length HILT from <i>sides</i> of backboard (Up to ten team members)			
Laid the backboard on top of a full length HILT.			
Deployed an extender for each team member that will lift the patient.			
Threaded each extender through the nearest backboard handhold.			
Packaged the patient according to local directives.			
Used correct lifting technique to minimize injury exposure to self.			
Lifted the patient in unison with team lift.			
Lowered the patient safely.			

Evaluator Name: _____

Date: _____