

Name: _____ Date: _____

Single Person Rescue



Single Person Rescue Practical

DIRECTIONS: Review the Single Person Rescue video in HILT U. Mark whether the candidate completed the task correctly or incorrectly.

Task	Correct	Incorrect	Notes
Folded one end of the HILT to the 2 nd set of extenders.			
Logrolled the patient, slid HILT into position, and logrolled patient onto HILT.			
Positioned the patient's head safely within the HILT.			
Deployed safety straps and fastened around the patient.			
Inserted hands into low lift base handles and dragged the patient while standing.			
Inserted hands into low lift base handles and dragged the patient while on knees.			

Evaluator Name: _____

Date: _____